


# Silicone Calibration Fluids for TOKI Viscometers


**Sintech Scientific (S.E.A.) Pte Ltd**  
 No. 194 Pandan Loop #04-11  
 Pantech Business Hub  
 Singapore 128383  
 Tel : 65-67790007  
 Fax : 63520877, 63520885  
 Email : [sales@sintechscientific.com](mailto:sales@sintechscientific.com)  
 Website : [www.sintechscientific.com](http://www.sintechscientific.com)

TOKI has also developed silicone calibration fluids which enable more efficient Calibration of our single cylinder rotational viscometer ( i.e. spindle rotor Viscometers such as our Model B ).



Type	Viscosity(mPa · s) 25°C	Main Application
<b>SC100</b>	90	M1 rotor BM II RB-85L TV-10M/15M TV-100BL/200BL
<b>SC2000</b>	1800	H2 rotor BH II RB-85H TV-10H/15H TV-100BH/200BH