

Laboratory Shaker (Horizontal Shaking)

MK201D


Rotary and reciprocate shaking motion

- Compact, space-saving design.
- Changeable rotary and reciprocate motion for mixing, extracting and stirring of samples.
- Digital display of shaking speed.
- Equipped with shaking timer.
- Easy to assemble and remove accessories.
- Various shaking modes when used with different shaking stages and racks.

Specifications

Product code	232024
Model	MK201D
Shaking mode	Rotary and Reciprocate (Manual operation)
Shaking width	Rotary: 30mm Reciprocate: 30mm
Shaking speed	20~200rpm, Stepless adjustment
Display type	Digital
Timer	0.5sec.~100hrs. (Unit and scale changeover type)
Allowable load weight	12kg
External dimensions*1	W442×D415×H130mm
Weight	19kg
Power source*2 (50/60Hz)	AC220V 1A Single phase with step-down transformer

*1 Protrusions not included. *2 The length of the power cord is about 2m outside the unit.



*Glassware is not included

Example of using a two-stage spring shaking rack (optional)

Operational accessories

Shaking rack for erlenmeyer flask clamp, Fixed type



Product code	Capacity	No. of clamps
232170	100mL	20pcs.
232171	200mL	10pcs.
232172	300mL	7pcs.
232173	500mL	4pcs.
232174	1,000mL	2pcs.

Shaking rack for erlenmeyer flask clamp, Inclined type



2 sets example (Inclination stage sold separately)

For 1 set		
Product code	Capacity	No. of clamps
232175	100mL	8pcs.
232176	200mL	4pcs.
232177	500mL	2pcs.

Shaking rack for spitz tube



2 sets example (Inclination stage sold separately)

For 1 set	
Product code: 232178	
Size 15ml	12pcs.
Applicable O.D. ø16.5×L110mm or less	

*Two sets can be shaken at the same time.

Shaking rack for centrifugal tube



2 sets example (Inclination stage sold separately)

For 1 set	
Product code: 232179	
Size 50mL	8pcs.
Applicable O.D. ø30×L110mm or less	

*Two sets can be shaken at the same time.

Shaking rack for SUS spring almighty type A



A type: for polyethylene bottle	
Product code: 232180	
Size 250ml	up to 9pcs.

Shaking rack for SUS spring almighty type B



B type: for erlenmeyer flask and test tube	
Product code: 232181	
17×15 = 255 grids	

Shaking rack with sticky non-skid sheet



Product code: 232182	
W400×D330mm	Material: silicon rubber
(For flat bottom flask use only)	

Two-stage shaking rack for SUS spring almighty type



C type: for erlenmeyer flask and test tube	
Product code: 231398	
For erlenmeyer flask	100ml×16pcs, 200mL×10pcs,
	500ml×4pcs, 1,000mL×2pcs.
For test tube ø16mm×105pcs.	when 45°inclination

Inclination setting stage for 232175~232179



Product code: 232078	
Necessary in tilt cases, up to 2 sets	

*Glassware and polyethylene bottles not included.