

Oil Bath B0302-E

[Page Top](#)[Product Summary](#)[Specifications](#)[Same series products](#)

Oil Bath

B0302-E

Item Code : 255426

Model : B0302-E

[Page Top](#)[Product Summary](#)[Specifications](#)[Same series products](#)

Oil Bath (B0302-E)

[Contact Us](#)[Printing](#)

Sintech Scientific (S.E.A.) Pte Ltd

No. 194 Pandan Loop #04-11

Pantech Business Hub

Singapore 128383

Tel : 65-67790007

Fax : 63520877, 63520885

Email : sales@sintechscientific.com

Website : www.sintechscientific.com

Product Summary

Features of Oil Bath (B0302-E)

Suitable for rotary evaporator RE202 series

Specifications

Specifications

Item Code	255426
Model	B0302-E
Operating temp. range	Room temp. +10 to 180°C

Oil BathB0302-E

Page Top	Product Summary	Specifications	Same series products
Temp. setting / display method	White LED digital display, key entry, minimum digit of 1oC		
Heater	1kW		
Sensor	Pt 100 ohm		
Safety device	Automatic overheat prevention, Independent overheat prevention (fixed temp.)		
Additional function	Overshoot alert, Auto resume (selectable), Calibration off-set		
External dimensions (φ×H mm)	φ262×240mm		
Remarks(dimensions)	Bath internal capacity : Approx.5L		
Power Source	AC100 to 115V 10 to 12A (50/60Hz) Service outlet 2A excluded		
Weight	4.5kg		

Maintenance and warranty

Conditions: temperature and humidity: 23 deg.C+/-5 deg.C, 65%RH+/-20% (no load)

External dimensions do not include protrusions.

Same series products



Economical Oil Bath (B0400)



Economical Oil Bath (B0410)



Economical Oil Bath (B0500)

