

Oil Bath B0302

[Page Top](#)[Product Summary](#)[Specifications](#)[Same series products](#)

Oil Bath B0302

Item Code : 255358

Model : B0302

[Page Top](#)[Product Summary](#)[Specifications](#)[Same series products](#)

Oil Bath (B0302)

[Contact Us](#)[Printing](#)**Sintech Scientific (S.E.A.) Pte Ltd**

No. 194 Pandan Loop #04-11

Pantech Business Hub

Singapore 128383

Tel : 65-67790007

Fax : 63520877, 63520885

Email : sales@sintechscientific.comWebsite : www.sintechscientific.com

Product Summary

Features of Oil Bath (B0302)

Suitable for rotary evaporator RE202 series

Specifications

Specifications

Item Code	255358
Model	B0302
Operating temp. range	Room temp. +10 to 180°C

Oil Bath B0302

[Page Top](#)[Product Summary](#)[Specifications](#)[Same series products](#)

Temp. setting / display method	White LED digital display, key entry, minimum digit of 1oC
Heater	1kW
Sensor	Pt 100 ohm
Safety device	Automatic overheat prevention, Independent overheat prevention (fixed temp.)
Additional function	Overshoot alert, Auto resume (selectable), Calibration off-set
External dimensions (φ×H mm)	φ262×240mm
Remarks(dimensions)	Bath internal capacity : Approx.5L
Power Source	AC100 to 115V 10A
Weight	4.5kg

Maintenance and warranty

Conditions: temperature and humidity: 23 deg.C+/-5 deg.C, 65%RH+/-20% (no load)

External dimensions do not include protrusions.

Same series products



Oil Bath (B0312-D)



Economical Oil Bath (B0400)



Economical Oil Bath (B0410)

