

Oil BathB0312-D

[Page Top](#)[Product Summary](#)[Specifications](#)[Same series products](#)

Oil Bath

B0312-D

Item Code : 255404

Model : B0312-D

[Page Top](#)[Product Summary](#)[Specifications](#)[Same series products](#)

Oil Bath (B0312-D)

[Contact Us](#)[Printing](#)

Sintech Scientific (S.E.A.) Pte Ltd

No. 194 Pandan Loop #04-11

Pantech Business Hub

Singapore 128383

Tel : 65-67790007

Fax : 63520877, 63520885

Email : sales@sintechscientific.com

Website : www.sintechscientific.com

Product Summary

Features of Oil Bath (B0312-D)

Suitable for rotary evaporator RE212 series

Specifications

Specifications

Item Code	255404
Model	B0312-D
Operating temp. range	Room temp. +10 to 180°C

Oil BathB0312-D

[Page Top](#)[Product Summary](#)[Specifications](#)[Same series products](#)

Temp. setting / display method	White LED digital display, key entry, minimum digit of 1oC
Heater	1kW
Sensor	Pt 100 ohm
Safety device	Automatic overheat prevention, Independent overheat prevention (fixed temp.)
Additional function	Overshoot alert, Auto resume (selectable), Calibration off-set
External dimensions (φ×H mm)	φ262×240mm
Remarks(dimensions)	Bath internal capacity : Approx.5L
Power Source	AC200 to 230V 5 to 6A (50/60Hz) Service outlet for RE212 : None
Weight	4.5kg

Maintenance and warranty

Conditions: temperature and humidity: 23 deg.C+/-5 deg.C, 65%RH+/-20% (no load)

External dimensions do not include protrusions.


Same series products



Economical Oil Bath (B0410)



Economical Oil Bath (B0500)



Economical Oil Bath (B0601)

